



Instrument Tip Protectors - 9.5 x 25.4mm/Orange



These tip protectors are a convenient, color-coded way to protect surgical instrument tips and edges. They provide exceptional protection and are vented to prevent rusting from condensation. The Assorted Pack (Item No. 29060-00) is 5 of each of the 6 available sizes.

General Info

Product Number	29060-06
Category	Instrument Care > Storage
EAN	4251909508500



Find 29060-06 on [finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	25 g	Net Weight	21 g
Package Size	25 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Plastisol		
Autoclavesafe	Yes		
Color	Orange		
Dimensionstext	9.5 x 25.4 mm		