



Instrument Tip Protectors - 15.9 x 25.4mm/Brown

These tip protectors are a convenient, color-coded way to protect surgical instrument tips and edges. They provide exceptional protection and are vented to prevent rusting from condensation. The Assorted Pack (Item No. 29060-00) is 5 of each of the 6 available sizes.



This product is discontinued and has limited availability.

General Info

Product Number	29060-07
Category	Instrument Care > Storage
EAN	4251909510015



Find 29060-07 on finescience.com

Technical Info

Shipping Weight	40 g	Net Weight	1 g
Package Size	25 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Plastisol		
Autoclavesafe	Yes		
Color	Brown		
Dimensionstext	15.9 x 25.4mm		