



Double Tip Protectors - Assorted Pack



These vented tip protectors are designed to hold hemostats and rongeurs firmly while allowing the jaws, box hinge and tips to stay open for thorough sterilization.

General Info

Product Number	29066-00
Category	Instrument Care > Storage
EAN	4251909508517



Find 29066-00 on [finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	15 g	Net Weight	8 g
Package Size	12 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Plastisol		
Autoclavesafe	Yes		
Color	Assorted		
Dimensionstext	Assorted		