



Double Tip Instrument Protectors - 12.7 x 20mm/Green

These vented tip protectors are designed to hold hemostats and rongeurs firmly while allowing the jaws, box hinge and tips to stay open for thorough sterilization.



General Info

Product Number	29066-01
Category	Instrument Care > Storage
EAN	4251909508524



Find 29066-01 on finescience.com

Technical Info

Shipping Weight	10 g	Net Weight	8 g
Package Size	10 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Plastisol		
Autoclavesafe	Yes		
Color	Green		
Dimensionstext	12.7 x 2 mm		