



Double Tip Instrument Protectors - 9 x 25.5mm/Pink

These vented tip protectors are designed to hold hemostats and rongeurs firmly while allowing the jaws, box hinge and tips to stay open for thorough sterilization.



General Info

|                |                           |
|----------------|---------------------------|
| Product Number | 29066-02                  |
| Category       | Instrument Care > Storage |
| EAN            | 4251909508531             |



Find 29066-02 on  
[finescience.com](https://finescience.com)

Technical Info

|                        |             |                 |     |
|------------------------|-------------|-----------------|-----|
| Shipping Weight        | 10 g        | Net Weight      | 8 g |
| Package Size           | 10 Pc.      | Unit of Measure | PA  |
| Minimum Order Quantity | 1           |                 |     |
| Alloymaterial          | Plastisol   |                 |     |
| Autoclavesafe          | Yes         |                 |     |
| Color                  | Pink        |                 |     |
| Dimensionstext         | 9 x 25.5 mm |                 |     |