

F · S · T

## Double Tip Instrument Protectors - 15 x 25mm/Yellow



These vented tip protectors are designed to hold hemostats and rongeurs firmly while allowing the jaws, box hinge and tips to stay open for thorough sterilization.

## General Info

Product Number	29066-03
Category	Instrument Care > Storage
EAN	4251909508548



Find 29066-03 on  
[finescience.com](https://www.finescience.com)

## Technical Info

Shipping Weight	20 g	Net Weight	18 g
Package Size	10 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloy material	Plastisol		
Autoclavesafe	Yes		
Color	Yellow		
Dimensionstext	15 x 25 mm		