



Double Tip Instrument Protectors - 25 x 26mm/Orange



These vented tip protectors are designed to hold hemostats and rongeurs firmly while allowing the jaws, box hinge and tips to stay open for thorough sterilization.

General Info

Product Number	29066-04
Category	Instrument Care > Storage
EAN	4251909508555



Find 29066-04 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	18 g
Package Size	10 Pc.	Unit of Measure	PA
Minimum Order Quantity	1		
Alloymaterial	Plastisol		
Autoclavesafe	Yes		
Color	Orange		
Dimensionstext	25 x 26 mm		