




Instrument Marking Tape - Blue/510cm

Organize instruments quickly and easily with this economical, non-toxic pressure sensitive tape. Color is unaffected by boiling, cold sterilization or autoclaving. 3 colors available. Autoclavable.



General Info

Product Number	29070-02
Category	Instrument Care > Marking
EAN	4251909510039



Find 29070-02 on finescience.com

Technical Info

Shipping Weight	10 g	Net Weight	6 g
Package Size	1 Roll	Unit of Measure	RO
Minimum Order Quantity	1		
Alloymaterial	Vinyl		
Autoclavesafe	Yes		
Color	Blue		
Lengthofroll	510 cm		
Tapewidth	0.6 cm		