




Instrument Marking Tape - Red/510cm

Organize instruments quickly and easily with this economical, non-toxic pressure sensitive tape. Color is unaffected by boiling, cold sterilization or autoclaving. 3 colors available. Autoclavable.



General Info

|                |                           |
|----------------|---------------------------|
| Product Number | 29070-04                  |
| Category       | Instrument Care > Marking |
| EAN            | 4251909510053             |



Find 29070-04 on [finescience.com](https://finescience.com)

Technical Info

|                        |        |                 |     |
|------------------------|--------|-----------------|-----|
| Shipping Weight        | 10 g   | Net Weight      | 6 g |
| Package Size           | 1 Roll | Unit of Measure | RO  |
| Minimum Order Quantity | 1      |                 |     |
| Alloymaterial          | Vinyl  |                 |     |
| Autoclavesafe          | Yes    |                 |     |
| Color                  | Red    |                 |     |
| Lengthofroll           | 510 cm |                 |     |
| Tapewidth              | 0.6 cm |                 |     |