



Instrument Marking Tape - White/510 cm

Tape Width: 0.6cm Length of Roll: 510cm Material: Vinyl Autoclave Safe: Yes Color: White Organize instruments quickly and easily with this economical, non-toxic pressure sensitive tape. Color is unaffected by boiling, cold sterilization or autoclaving. 9



General Info

Product Number	29070-05
Category	Instrument Care > Marking
EAN	4251909512354



Find 29070-05 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	10 g	Net Weight	6 g
Package Size	1 Roll	Unit of Measure	RO
Minimum Order Quantity	1		
Alloymaterial	Vinyl		
Autoclavesafe	Yes		
Color	White		
Lengthofroll	510 cm		
Tapewidth	0.6 cm		