



Instrument Marking Tape - Brown/290cm

Organize instruments quickly and easily with this economical, non-toxic, pressure sensitive tape. Color is unaffected by boiling, cold sterilization or autoclaving. 7 colors available. Autoclavable.



This product is discontinued and has limited availability.

General Info

Product Number	29071-07
Category	Instrument Care > Marking
EAN	4251909508623



Find 29071-07 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	6 g
Package Size	1 Roll	Unit of Measure	RO
Minimum Order Quantity	1		
Alloymaterial	Vinyl		
Autoclavesafe	Yes		
Color	Brown		
Lengthofroll	290 cm		
Tapewidth	0.6 cm		