



Instrument Marking Tape - Orange/290cm

Organize instruments quickly and easily with this economical, non-toxic, pressure sensitive tape. Color is unaffected by boiling, cold sterilization or autoclaving. 7 colors available. Autoclavable.



General Info

Product Number	29071-08
Category	Instrument Care > Marking
EAN	4251909508630



Find 29071-08 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	20 g	Net Weight	6 g
Package Size	1 Roll	Unit of Measure	RO
Minimum Order Quantity	1		
Alloymaterial	Vinyl		
Autoclavesafe	Yes		
Color	Orange		
Lengthofroll	290 cm		
Tapewidth	0.6 cm		