



Instrument Marking Tape - 7 assorted colors/107cm each



Organize instruments quickly and easily with this economical, non-toxic, pressure sensitive tape. Color is unaffected by boiling, cold sterilization or autoclaving. 7 assorted colors: Yellow, Red, Black, White, Green, Blue and Orange.

General Info

Product Number	29080-02
Category	Instrument Care > Marking
EAN	4251909515072



Find 29080-02 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	17 g
Package Size	7 Rolls	Unit of Measure	RO
Minimum Order Quantity	1		
Alloymaterial	Vinyl		
Autoclavesafe	Yes		
Color	Assorted		
Lengthofroll	107 cm		
Tapewidth	0.6 cm		