



Wide screen LED display with the ability to measure in millimeters, inches, and fractions. Measure inside, outside, step and depth. Battery, spare battery and plastic storage case included.

General Info

Product Number	30087-00
Category	Lab Accessories > Measurement
EAN	4251909508692



Find 30087-00 on
finescience.com

Technical Info

Shipping Weight	300 g	Net Weight	295 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	23.5 cm		
Increments	0.1 mm		
Maximummeasurement	15.2 cm / 6 inch		