



Tough polymer is unbreakable under normal conditions. Measures inside, outside, step and depth, metric and inch to 0.1mm or 1/64. Capacity is 13cm.

General Info

Product Number	30087-20
Category	Lab Accessories > Measurement
EAN	4251909508708



Find 30087-20 on
finescience.com

Technical Info

Shipping Weight	105 g	Net Weight	42 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Plastic		
Autoclavesafe	Yes		
Length	21 cm		
Increments	0.1 mm		
Maximummeasurement	19 cm / 7 inch		