



The scissor blades are made of high quality stainless-steel and guarantee a perfect cut. Handles are made of glass fiber reinforced polyamide. Not autoclavable.

General Info

Product Number	37500-05
Category	Scissors & Micro Scissors > Utility Scissors
EAN	4251909508739



Find 37500-05 on
finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	21 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Stainless Steel (Inox) with plastic handle		
Autoclavesafe	No		
Cuttingedge	50 mm		
Length	13 cm		
Tips	Sharp/Sharp		
Tipshape	Straight		