



Utility Scissors - Straight/ Blunt-Blunt / 21 cm



These scissors are extremely hard 54 on the Rockwell scale with the ability to cut tin, aluminum, plastic, bone and other hard materials. The strong 3mm thick blades are magnetic, and the handles are plastic. Not autoclavable.

This product is discontinued and has limited availability.

General Info

Product Number	37501-00
Category	Scissors & Micro Scissors > Utility Scissors
EAN	4251909511470

Technical Info

Shipping Weight	100 g	Net Weight	82 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Stainless Steel (Inox) with plastic handle		
Cuttingedge	80 mm		
Length	21 cm		
Sterile	No		
Tips	Blunt/Blunt		
Tipshape	Straight		