



Vascular Catheterization Surgical Pack - For models between 100g and 400g

This surgical kit is designed for performing vascular catheterization procedures on models between 100g and 400g commonly used in biomedical research. It contains all the necessary instruments to conduct precise and delicate surgeries with ease especially implanting Vascular Access Buttons™ by Instech Laboratories. Fine Science Tools highly recommends Instech's revolutionary Vascular Access Buttons™ featuring PinPort™ technology that simply and cleanly connects an implanted catheter to a syringe for direct access, or to a tether for continuous access. An industry leading refinement over many other methods of vascular access, Instech Laboratories is quickly becoming the preferred method of catheter exteriorization.



Content of kit:

- 1 x Vannas-Tübingen Spring Scissors - Angled Up/Sharp/9.5cm/5mm Cutting Edge [15005-08]
- 1 x Strabismus Scissors - Straight/Blunt-Blunt/11.5cm [14074-11]
- 1 x Wagner Scissors - Straight/Sharp-Sharp/12cm [14068-12]
- 2 x S&T Vessel Dilating Forceps- Angled 45°/13.5cm/0.2mm Tip Diameter [00276-13]
- 1 x Dumont 545 Forceps - Micro-Blunted Standard Tips/Angled 45°/Inox/11cm [11253-25]
- 1 x Graefe Forceps - 1x2 Teeth/Curved/10cm [11054-10]
- 1 x Kalt Needle Holder, straight, 13.5cm, cross-serrated and serrated [12080-14]
- 1 x Micro-Mosquito Hemostat - Serrated/Straight/12cm [13010-12]
- 1 x Grooved Director with Tongue, 13 cm [10290-13]
- 1 x Schwartz Micro Serrefine - Light Bend/26mm [18052-02]
- 1 x Magnetic Fixator Retraction System - Elastomer [18200-07]
- 1 x Magnetic Fixator Retraction System - Blunt Retractors/5mm Wide [18200-11]
- 1 x Absorption Spears - Mounted/Sterile [18105-01]
- 1 x Instrument Tip Protectors - Assorted Pack [29060-00]
- 1 x Suture Thread - Braided Silk/Size 6/0 [18020-60]
- 1 x Plastic Sterilization Container-Single Tier - 25 x 15 x 2cm [20840-03]

General Info

Product Number	70003-02
Category	Kits
EAN	4251909511715



Find 70003-02 on
finescience.com

Technical Info

Shipping Weight	1010 g	Net Weight	802 g
Package Size	1 Pc.	Unit of Measure	SET
Minimum Order Quantity	1		
Feature	For models between 100g and 400g		