



Professional Zebrafish Kit



The Professional Zebrafish Kit is ideal for advanced users and contains all the essential instruments for the most common procedures. It comes in an autoclavable box that can also be used for storage.

- Content of kit:
- 1 x Probe - Straight/15.5cm/0.25mm Tip Diameter [10140-01]
 - 1 x Probe - Angled 80° Short/15.5cm/0.15mm Tip Diameter [10140-02]
 - 1 x Tissue Forceps - 1x2 Teeth/Straight/12cm [11021-12]
 - 1 x Graefe Forceps - Serrated/Curved/10cm [11051-10]
 - 1 x Dumont #5 Forceps - Standard Tips/Straight/Inox/11cm [11251-20]
 - 1 x Noyes Spring Scissors - Straight/Sharp/12cm/14mm Cutting Edge [15012-12]
 - 1 x Plastic Sterilization Container-Single Tier - 17 x 10 x 2cm [20840-02]

General Info

Product Number	70004-01
Category	Kits
EAN	4251909511487



Find 70004-01 on
finescience.com

Technical Info

Shipping Weight	550 g	Net Weight	540 g
Package Size	1 Pc.	Unit of Measure	SET
Minimum Order Quantity	1		