



Zebrafish Student Kit



- Inhalt des Kits:
- 1 x Probe - Straight/15.5cm/0.25mm Tip Diameter [10140-01]
 - 1 x Probe - Angled 80° Short/15.5cm/0.15mm Tip Diameter [10140-02]
 - 1 x Graefe Forceps - Serrated/Angled/10cm [11049-10]
 - 1 x Dumont #5 Forceps - Standard Tips/Straight/Inox/11cm [11251-20]
 - 1 x Student Vannas Spring Scissors - Straight/Sharp/9cm/5mm Cutting Edge [91500-09]
 - 1 x Plastic Sterilization Container-Single Tier - 17 x 10 x 2cm [20840-02]

FST Student Instruments wurden mit Blick auf Qualität und Wirtschaftlichkeit entwickelt, daher verfügen sie über weniger manuelle Nachbearbeitung und werden aus kostengünstigeren Legierungen hergestellt.

Allgemein

Produktnummer	70004-02
Kategorie	Student Line > Kits
EAN	4251909511647



Finde 70004-02 auf finescience.com

Technische Informationen

Versandgewicht	530 g	Nettogewicht	520 g
Paketgröße	1 Pc.	Maßeinheit	SET
Mindestbestellmenge	1		
Feature	Student Quality		