



## Student Anatomical Standard Forceps - Straight/Student Stainless Steel/13cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

### General Info

Product Number	91100-13
Category	Student Line > Forceps & Tweezers
EAN	4251909510947



Find 91100-13 on  
[finescience.com](https://www.finescience.com)

### Technical Info

Shipping Weight	20.5 g	Net Weight	19.5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloy material	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	13 cm		
Feature	Student Quality		
Tipwidth	2.7 mm		
Tips	Serrated		
Tipshape	Straight		