



Student Semken Forceps - Serrated/Straight/Student Stainles Steel/11.5cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91108-12
Category	Student Line > Forceps & Tweezers
EAN	4251909511012



Find 91108-12 on finescience.com

Technical Info

Shipping Weight	16.5 g	Net Weight	15.5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11.5 cm		
Feature	Student Quality		
Tipwidth	1.2 mm		
Tips	Serrated		
Tipshape	Straight		