



Student #5/45 - 45° Angled/Student Stainless Steel/11cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91110-10
Category	Student Line > Forceps & Tweezers
EAN	4251909508746



Find 91110-10 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	13 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11 cm		
Feature	Student Quality		
Style	#5/45		
Tipdimensions	0.2 x 0.2 mm		
Tips	Standard		
Tipshape	Angled		