



Student Micro-Adson Forceps - 1x2 Teeth/Straight/Student Stainless Steel/12cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91119-12
Category	Student Line > Forceps & Tweezers
EAN	4251909511159



Find 91119-12 on finescience.com

Technical Info

Shipping Weight	24.5 g	Net Weight	23.5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	12 cm		
Feature	Student Quality		
Tipwidth	1.3 mm		
Tips	1x2 Teeth		
Tipshape	Straight		