



Student Surgical Standard Forceps - Straight/Student Stainless Steel/12cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91121-12
Category	Student Line > Forceps & Tweezers
EAN	4251678878422



Find 91121-12 on finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	14 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	12 cm		
Feature	Student Quality		
Tipwidth	2.2 mm		
Tips	1x2 Teeth		
Tipshape	Straight		