



## Student Surgical Standard Forceps - Straight/Student Stainless Steel/14.5cm

Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.



### General Info

Product Number	91121-14
Category	Student Line > Forceps & Tweezers
EAN	4251909511067



Find 91121-14 on  
[finescience.com](https://www.finescience.com)

### Technical Info

Shipping Weight	24 g	Net Weight	23 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloy material	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	14.5 cm		
Feature	Student Quality		
Tipwidth	2.6 mm		
Tips	1x2 Teeth		
Tipshape	Straight		