



Student Surgical Standard Forceps - Straight/Student Stainless Steel/20cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91121-20
Category	Student Line > Forceps & Tweezers
EAN	4251909511104



Find 91121-20 on finescience.com

Technical Info

Shipping Weight	46 g	Net Weight	45 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	20 cm		
Feature	Student Quality		
Tipwidth	4 mm		
Tips	1x2 Teeth		
Tipshape	Straight		