



Student Surgical Narrow Forceps - Straight/Student Stainless Steel/16cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91123-16
Category	Student Line > Forceps & Tweezers
EAN	4251909511142



Find 91123-16 on finescience.com

Technical Info

Shipping Weight	27.5 g	Net Weight	26.5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	16 cm		
Feature	Student Quality		
Tipwidth	2.4 mm		
Tips	1x2 Teeth		
Tipshape	Straight		