



Student Gillies Forceps - 1x2 Teeth & Cross Serrations/Straight/Student
Stainless Steel/15.5cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91128-15
Category	Student Line > Forceps & Tweezers
EAN	4251909511180



Find 91128-15 on
finescience.com

Technical Info

Shipping Weight	28 g	Net Weight	27 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	15.5 cm		
Feature	Student Quality		
Tipwidth	2.1 mm		
Tips	1x2 Teeth and Cross-Serrations		
Tipshape	Straight		