



Student Dumont #5 Forceps - Standard Tips/Straight/Inox/11cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91150-20
Category	Student Line > Forceps & Tweezers Most Popular > Student Line
EAN	7640241754268



Find 91150-20 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	11 cm		
Feature	Student Quality		
Style	#5		
Tipdimensions	0.1 x 0.06 mm		
Tips	Standard		
Tipshape	Straight		