



Student Splinter Forceps - Serrated/Straight/Student Stainless Steel/10.5cm/0.3mm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91156-11
Category	Student Line > Forceps & Tweezers
EAN	4251909511203



Find 91156-11 on finescience.com

Technical Info

Shipping Weight	30 g	Net Weight	10 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	10.5 cm		
Feature	Student Quality		
Tipdiameter	0.3 mm		
Tipwidth	0.3 mm		
Tips	Serrated		
Tipshape	Straight		