



Student Allis Clamp - 5x6 Teeth/Straight/Student Stainless Steel/15cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91191-15
Category	Student Line > Clamps & Hemostats
EAN	4251909511272



Find 91191-15 on
finescience.com

Technical Info

Shipping Weight	26.5 g	Net Weight	25.5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	15 cm		
Feature	Student Quality		
Tipwidth	6 mm		
Tips	5x6 Teeth		
Tipshape	Straight		