



Student Tohoku Towel Clamp - Serrated/Student Stainless Steel/10cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91194-10
Category	Student Line > Clamps & Hemostats
EAN	4251909510893



Find 91194-10 on
finescience.com

Technical Info

Shipping Weight	24.5 g	Net Weight	23.5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	10 cm		
Feature	Student Quality		
Tipwidth	2.5 mm		
Tips	Round and serrated tips		
Tipshape	Curved		