



Student Hose Clamp - Cross Serrated/Straight/Student Stainless Steel/15cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91198-15
Category	Student Line > Clamps & Hemostats
EAN	4251909510923



Find 91198-15 on finescience.com

Technical Info

Shipping Weight	42 g	Net Weight	41 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	15 cm		
Feature	Student Quality		
Tipwidth	4 mm		
Tips	Cross-Serrated		
Tipshape	Straight		