



Student Mayo-Hegar Needle Holder - Serrated/Straight/Student Stainless Steel/16cm/with Lock



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

|                |                               |
|----------------|-------------------------------|
| Product Number | 91204-16                      |
| Category       | Student Line > Needle Holders |
| EAN            | 4251909511265                 |



Find 91204-16 on [finescience.com](https://finescience.com)

Technical Info

|                        |                                 |                 |      |
|------------------------|---------------------------------|-----------------|------|
| Shipping Weight        | 41 g                            | Net Weight      | 40 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA   |
| Minimum Order Quantity | 1                               |                 |      |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |      |
| Clampinglength         | 21 mm                           |                 |      |
| Length                 | 16 cm                           |                 |      |
| Feature                | Student Quality                 |                 |      |
| Lockable               | Yes                             |                 |      |
| Scissors               | No                              |                 |      |
| Serrated               | Yes                             |                 |      |
| Tipwidth               | 2.6 mm                          |                 |      |
| Tipshape               | Straight                        |                 |      |