



Student Crile-Rankin Hemostat - Serrated/Curved/Student Stainless Steel/16cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

|                |                                   |
|----------------|-----------------------------------|
| Product Number | 91305-16                          |
| Category       | Student Line > Clamps & Hemostats |
| EAN            | 4251909510824                     |



Find 91305-16 on [finescience.com](https://finescience.com)

Technical Info

|                        |                                 |                 |        |
|------------------------|---------------------------------|-----------------|--------|
| Shipping Weight        | 38.5 g                          | Net Weight      | 37.5 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA     |
| Minimum Order Quantity | 1                               |                 |        |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |        |
| Length                 | 16 cm                           |                 |        |
| Feature                | Student Quality                 |                 |        |
| Tipwidth               | 2 mm                            |                 |        |
| Tips                   | Serrated                        |                 |        |
| Tipshape               | Curved                          |                 |        |