



Student Halsted-Mosquito Hemostat - Serrated/Straight/Student
Stainless Steel/12.5cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91308-12
Category	Student Line > Clamps & Hemostats
EAN	4251678878545



Find 91308-12 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	21 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Clampinglength	20 mm		
Length	12.5 cm		
Feature	Student Quality		
Tipwidth	1.4 mm		
Tips	Serrated		
Tipshape	Straight		