



Student Rochester-Pean Hemostat - Serrated/Straight/Student Stainless Steel/16cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91316-16
Category	Student Line > Clamps & Hemostats
EAN	4251909510664



Find 91316-16 on finescience.com

Technical Info

Shipping Weight	44.5 g	Net Weight	43.5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Length	16 cm		
Feature	Student Quality		
Tipwidth	2.5 mm		
Tips	Serrated		
Tipshape	Straight		