



Student Kelly Hemostat - Serrated/Curved/Student Stainless Steel/14cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91319-14
Category	Student Line > Clamps & Hemostats
EAN	4251909510787



Find 91319-14 on
finescience.com

Technical Info

Shipping Weight	30.5 g	Net Weight	29.5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	14 cm		
Feature	Student Quality		
Tipwidth	2 mm		
Tips	Serrated		
Tipshape	Curved		