



Student Halsted-Mosquito Hemostat - 1x2 Teeth and Serrations/Curved/Student Stainless Steel/12.5cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91331-12
Category	Student Line > Clamps & Hemostats
EAN	4251909510831



Find 91331-12 on finescience.com

Technical Info

Shipping Weight	21 g	Net Weight	20 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Length	12.5 cm		
Feature	Student Quality		
Tipwidth	2 mm		
Tips	1x2 Teeth and Serrations		
Tipshape	Curved		