



Student Standard Surgical Scissors - Straight/Sharp-Blunt/Student
Stainless Steel/10.5cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91401-10
Category	Student Line > Scissors & Micro Scissors
EAN	4251678877968



Find 91401-10 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	21 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	30 mm		
Cuttingedge	20 mm		
Length	10.5 cm		
Feature	Student Quality		
Serrated	No		
Tips	Sharp/Blunt		
Tipshape	Straight		