



Student Standard Surgical Scissors - Straight/Sharp-Sharp/Student
Stainless Steel/13cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91402-13
Category	Student Line > Scissors & Micro Scissors
EAN	4251909510442



Find 91402-13 on
finescience.com

Technical Info

Shipping Weight	36.5 g	Net Weight	35.5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	42 mm		
Cuttingedge	28 mm		
Length	13 cm		
Feature	Student Quality		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Straight		