



Student Standard Surgical Scissors - Curved/Sharp-Blunt/Student
Stainless Steel/12cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91404-12
Category	Student Line > Scissors & Micro Scissors
EAN	4251909510459



Find 91404-12 on
finescience.com

Technical Info

Shipping Weight	28.5 g	Net Weight	27.5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	39 mm		
Cuttingedge	26 mm		
Length	12 cm		
Feature	Student Quality		
Serrated	No		
Tips	Sharp/Blunt		
Tipshape	Curved		