



Student Standard Surgical Scissors - Curved/Sharp-Sharp/Student Stainless Steel/12cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91405-11
Category	Student Line > Scissors & Micro Scissors
EAN	4251909510480



Find 91405-11 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	31 g	Net Weight	30 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Bladlength	39 mm		
Cuttingedge	26 mm		
Length	12 cm		
Feature	Student Quality		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Curved		