



Student Standard Surgical Scissors - Curved/Sharp-Sharp/Student
Stainless Steel/14.5cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91405-14
Category	Student Line > Scissors & Micro Scissors
EAN	4251909510503



Find 91405-14 on
finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	49 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Bladelength	45 mm		
Cuttingedge	30 mm		
Length	14.5 cm		
Feature	Student Quality		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Curved		