



Student Mayo Scissors - Curved/Blunt-Blunt/Student Stainless Steel/17cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91411-17
Category	Student Line > Scissors & Micro Scissors
EAN	4251909510534



Find 91411-17 on finescience.com

Technical Info

Shipping Weight	68.5 g	Net Weight	67.5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	48 mm		
Cuttingedge	32 mm		
Length	17 cm		
Feature	Student Quality		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Curved		