



Student Lister Scissors - Angled/Blunt-Blunt/Student Stainless Steel/14cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91413-14
Category	Student Line > Scissors & Micro Scissors
EAN	4061706159627



Find 91413-14 on finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	46 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	39 mm		
Cuttingedge	26 mm		
Length	14 cm		
Feature	Student Quality		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Angled to side		