



Student Lister Scissors - Angled/Blunt-Blunt/Student Stainless Steel/18cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91413-18
Category	Student Line > Scissors & Micro Scissors
EAN	4251909510626



Find 91413-18 on finescience.com

Technical Info

Shipping Weight	89.5 g	Net Weight	88.5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Bladelength	45 mm		
Cuttingedge	30 mm		
Length	18 cm		
Feature	Student Quality		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Angled to side		