



Student Metzenbaum Baby Scissors - Curved/Blunt-Blunt/Student
Stainless Steel/14.5cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

Product Number	91419-14
Category	Student Line > Scissors & Micro Scissors
EAN	4251909510541



Find 91419-14 on
finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	34 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	33 mm		
Cuttingedge	22 mm		
Length	14.5 cm		
Feature	Student Quality		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Curved		