



Student Adson Forceps - 1x2 Teeth & Cross Serrations/Straight/Student  
Stainless Steel/12cm



Our student instruments are designed with both quality and affordability in mind, featuring less hand-finishing and manufactured with less expensive alloys.

General Info

|                |                                   |
|----------------|-----------------------------------|
| Product Number | 91427-12                          |
| Category       | Student Line > Forceps & Tweezers |
| EAN            | 4251909511166                     |



Find 91427-12 on  
[finescience.com](https://finescience.com)

Technical Info

|                        |                                 |                 |      |
|------------------------|---------------------------------|-----------------|------|
| Shipping Weight        | 21 g                            | Net Weight      | 20 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA   |
| Minimum Order Quantity | 1                               |                 |      |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |      |
| Length                 | 12 cm                           |                 |      |
| Feature                | Student Quality                 |                 |      |
| Tipwidth               | 1.3 mm                          |                 |      |
| Tips                   | 1x2 Teeth/Cross-serrations      |                 |      |
| Tipshape               | Straight                        |                 |      |